


# **Mason Tayler Medical**

## **Gastric Lavage Tubes**

**. 100% Latex Free**

**. Cost Effective**

**. Radiopaque**



**Mason Tayler Medical**  
**888-270-7494**

# Gastric Lavage Tubes



- **Latex Free**
- **More efficient irrigation-multiple eyes and large smooth open distal tip.**
- **Easy intubation-satin finish surface provides less surface contact.**
- **Markings at 18", 22", 26" and 30".**
- **Clean ready to use, disposable non-sterile.**
- **Length of 48"**
- **Radiopaque.**

Item #	Size	Units
MTM 528	28 Fr.	10/box
MTM 532	32 Fr.	10/box
MTM 536	36 Fr.	10/box

Visit our Website at:  
[www.masontayler.com](http://www.masontayler.com)

For additional information or to place an order,  
please feel free to contact:

1-888-270-7494 / 716-883-8097  
or Fax: 716-883-8099

Email: [masontayler@worldnet.att.net](mailto:masontayler@worldnet.att.net)



MEDICAL PRODUCTS CORP.  
1083 Delaware Ave., Buffalo, NY 14209 USA



We accept  
Visa and Mastercard

