


Infu-Stat™

Disposable Pressure Infuser

- **100% Latex Free**
- **Cost Competitive**
- **User Friendly**
- **One hand operation**




Mason Tayer Medical
888-270-7494

Infu-Stat™

Disposable Pressure Infuser



- Latex Free
- Loads easily.
- One hand operation for inflation or deflation.
- Large bore stopcock for rapid inflation and deflation.
- Color-coded gauge is easy to read and maintains reliable pressure.
- The safety valve eliminates the possibility of over inflation and releases air at 325mm/hg.
- Sizing is color-coded for easy identification-Graphics are sharp and easy to read.
- Economical, single patient use to prevent cross contamination and eliminates costly cleaning and repair.



Item #	Size	Units
MTM 305	500ml	5/box - 25/cs
MTM 310	1000ml	5/box - 25/cs
MTM 330E	3000ml	2/box - 10/cs

Visit our Website at:
www.masontaylor.com

For additional information or to place an order,
please feel free to contact:

1-888-270-7494 / 716-883-8097
or Fax: 716-883-8099

Email: masontaylor@worldnet.att.net

We Gladly Accept:



MEDICAL PRODUCTS CORP.
1083 Delaware Ave., Buffalo, NY 14209 USA

