




LimbSecur

Body Positioners

. 100% Latex Free

. Cost Effective

. Reusable




Mason Taylor Medical


888-270-7494

LimbSecur

Body Positioners



MTE-902
Head Cradle



MTE-907-7" Donut
MTE-909-9" Donut

- **Latex Free**
- Long term reusability and cost effectiveness.
- **Will not block X-Rays.**
- **Moisture resistant, vapor-permeable, non-skid outer skin.**
- **Exclusive polyurethane foam that provides resilience use after use.**
- **Reduce disposal costs.**
- **Non-toxic, hypoallergenic and non-irritating material.**
- **Maximum economy and minimum waste.**
- **Clean with soap and water or cold sterilants.**
- **Can be flash autoclaved according to normal hospital protocols.**
- **Packaged individually.**

*For additional information or to place an order,
please feel free to contact:*

*Visit our Website at:
www.masontayler.com*



MEDICAL PRODUCTS CORP.
1083 Delaware Ave., Buffalo, NY 14209 USA

Mason Tayler Medical
1-888-270-7494 / 716-883-8097
or Fax: 716-883-8099
Email: masontayler@worldnet.att.net



We accept
Visa and Mastercard

