

# World's First Patented Antimicrobial Stethoscope

## Protect the lives you touch.

The new **MicrobeCare™** Stethoscopes are constructed with the **highest** quality materials and are the only **fully protected** antimicrobial stethoscopes available today, including the diaphragm, head, tubing, ear buds and all metal surfaces!



for better care...  
**MicrobeCare™**



**89%** Bottom Line Health reported that out of 150 physician stethoscopes tested, **133 were found to harbor staph bacteria.**

Read More





# Antimicrobial Stethoscope


for better care...

**MicrobeCare**

**MicrobeCare™** stethoscopes are **fully antimicrobial from top to bottom**.  
Manufactured with *superior materials* for improved acoustic performance  
and *long lasting durability*.

## Eartips

**Soft-silicon threaded eartips** compress to fit the ear canal  
for a superior acoustic seal and **unmatched comfort**.  
We include an extra pair with each scope.



## Chestpiece

The chestpiece is the portion of the stethoscope  
that conducts sound. The flat side is known as the diaphragm and  
is used to conduct sounds when placed on the patient.  
The chestpiece may have one diaphragm (double/single) or two diaphragms.

## Tubing

Made from latex-free PVC, our cardiology tubing has two distinct sound  
channels for improved stereo separation. Our dense, co-extruded material  
improves acoustic transmission by reducing ambient noise. Available in a  
number of **popular colors** for personalization!

## Auscultation Quality

Sound as good as it looks. **Superior materials**  
make for top quality precision acoustics.



The **MicrobeCare™** stethoscope lets you  
hear what you've been missing.



# How does it **really** work though?

At the molecular level (very small), **MicrobeCare™** works in simple terms by acting like a type of **sword** or **lance**. When microbes come into contact with the molecular coating on our stethoscopes, their cell-wall is **breached** by our patented **MicrobeCare™** technology. Thus, killing the microbe without degrading the quality of the antimicrobial coating for long lasting protection.

**MicrobeCare™** attracts the microbe to the lance.



conducts  
e membrane  
biology models  
ams (double/double).



The lance bursts the membrane, **killing the microbe.**

(without degrading the anti-microbial)

Ordering Information 

# Order Information

For color selections please visit: [www.Antimicrobialstethoscope.com/photos/](http://www.Antimicrobialstethoscope.com/photos/)

## Cardiology

All cardiology scopes ship in a large box



(CDS = Cardiology, Double Bell, Single Diaphragm)

Select: | CDS-BL    CDS-P    CDS-O    CDS-G    CDS-B    CDS-M

(CDD = Cardiology, Double Bell, Double Diaphragm)

Elite: | CDD-BL    CDD-P    CDD-O    CDD-G    CDD-B    CDD-M

(CSS = Cardiology, Single Bell, Single Diaphragm)

Supreme: | CSS-BL    CSS-P    CSS-O    CSS-G    CSS-B    CSS-M

(BL=Black, P=Pink, O=Orange, G=Green, B=Blue, M=Magenta)

## Legacy

All legacy scopes ship in a small box



(LDS = Legacy, Double Bell, Single Diaphragm)

Focus: | LDS-BL    LDS-P    LDS-O    LDS-G    LDS-B    LDS-M

(LPDS = Legacy Pediatric/Infant, Double Bell, Single Diaphragm)

Pediatric: | LPDS-BL    LPDS-P    LPDS-O    LPDS-G    LPDS-B    LPDS-M

(BL=Black, P=Pink, O=Orange, G=Green, B=Blue, M=Magenta)

\*All stethoscopes include accessory kit with additional diaphragm and earbuds.



1602 Barclay Blvd  
Bualo Grove, IL 60089  
[www.ParasolMed.com](http://www.ParasolMed.com)  
(847) 520-3000

FDA Registered Manufacturer  
FDA Listed Products  
ETL Listed UL Approved  
EMC Tested

Distributed By: