


CLARISONUS PLUS™

Disposable Esophageal Stethoscope with Listening Device



Clarisonus Plus is a new patent pending disposable Esophageal Stethoscope coupled with a disposable aural monoscope. Standard Esophageal Stethoscopes require the anesthesia care provider to attach an entirely separate disposable or reusable device to listen to the patient's heart and lung sounds during anesthesia. With the Clarisonus Plus, all you need to monitor the patient's core body temperature and listen to heart and lung sounds is right at your fingertips, from package to placement the product is all in one, making it simple to comply with Basic Anesthesia Monitoring Standards. Clarisonus Plus is a simple to use, cost effective approach to anesthesia monitoring.

ORDERING INFORMATION

Box Quantity : 30 each per box
Product is sold clean, non-sterile.

Catalog Number	Description	Size
NS4009-ESC	Clarisonus Plus Esophageal Stethoscope	9 Fr
NS40012-ESC	Clarisonus Plus Esophageal Stethoscope	12 Fr
NS40018-ESC	Clarisonus Plus Esophageal Stethoscope	18 Fr
NS40024-ESC	Clarisonus Plus Esophageal Stethoscope	24 Fr



22805 - B Savi Ranch Pkwy
Yorba Linda, CA 92887

PHONE: (714) 283-3099

FAX: (714) 283-3033

WEB SITE: www.starboardmedical.com

EMAIL: info@starboardmedical.com