



PARASOL

Patient Fall Prevention System

The only fall prevention system to address incontinence
the number one reason for falls

drysmart® incontinence pad real-time incontinence monitoring

The new drysmart® pad is part of a total fall prevention system focusing on elimination related falls. Real-time incontinence monitoring engineered into our unique, highly absorbent drysmart® pads notify clinicians when a patient has evacuated in bed. This early warning can allow intervention before the patient attempts to exit the bed.




- Collects up to 700 ml of fluid
- Superior wicking material helps draw moisture away from patient
- Easy connection & disconnection from the Parasol Medical™ Monitors
- Available for both wired & wireless Parasol Medical™ Monitors
- Integrates with existing nurse call systems

**“Urinary incontinence is associated
with an increase in falls”¹**

50%

According to peer reviewed healthcare studies elimination related falls are responsible for more than 50% of falls within acute care facilities.





combo bed/chair pad

The Parasol Medical combo pad combines pressure sensitive fall prevention & real-time incontinence monitoring, integrated into one convenient pad.



potty pad

The Parasol Medical potty pad is a key tool to the Parasol Fall Prevention System. As part of a total rooms solution clinicians can monitor patients using the commode or toilet.



led lights

Parasol has engineered its new wireless falls prevention monitor with strong LED lights to help patients navigate the path from bed to bathroom.

“The frequency of night time journeys to the toilet, often combined with poor lighting can increase the odds of falling significantly.”¹

Incontinent Pads & Accessories

- drySmart® Incontinence Pad WDS-01
- Wired Potty Pad Sensor 30 Day SP6S3
- Wireless Potty Pad Sensor 30 Day WP30
- Bed Pad Sensor Fall / Incontinence Combo.... IFCC8
- Chair Pad Sensor Fall / Incontinence Combo .. IFCC8
- drySmart® Wired Incontinence Quick Connector QC-01
- drySmart® Wireless Incontinence Quick Connector QC-02



1602 Barclay Blvd
Buffalo Grove, IL 60089
www.parasolmed.com
(847) 520-3000

FDA Registered Manufacturer
FDA Listed Products
ETL Listed UL Approved
EMC Tested



Distributed By: