

## Parasol Advanced Wireless Fall Prevention System

The first completely wireless, total-room solution that combines fall monitoring with real-time incontinence monitoring.



**PARASOL**

Treated with MicroBECARE™

# Prevent falls from every surface patients occupy with Parasol Medical's total room solution.



The Parasol Medical Wireless Patient Fall Prevention System is the first wireless patient fall prevention aid that combines fall monitoring with real-time incontinence monitoring, for proactive fall prevention.

This total room solution provides caregivers with the most advanced tools to prevent costly injuries due to falls and incontinence. Our patented wireless

technology eliminates cords and tripping hazards, and lets you configure patient rooms without restrictions due to cord length or proximity to nurse call outlets.

With its many powerful features—including status on demand that confirms connected and patient activated devices—it's the most advanced way to reduce your organization's exposure to risk and protect your patients like never before.



The Parasol Medical Wireless Patient Fall Prevention System includes a Wireless Bed Pad, Wireless Chair Pad, Wireless Potty or Commode Pad, Wireless Mobility Monitor with Voice Recorder (shown with wall bracket), Monitor Key, Nurse Call Interface, drysmart® Wireless Incontinence Pads, drysmart™ Quick Connector and a Wireless Seat Belt for wheelchair use.

## Revolutionary, Real-Time Incontinence Monitoring

Elimination-related falls are responsible for more than **50%** of falls within acute care facilities.\*

drysmart® incontinence pads monitor wetness in real-time and give nurses more time to reach patients before they get up and risk a fall. Each highly-absorbent pad is coated with moisture-sensitive thermal conductive ink, and sends an alarm in real-time when a patient evacuates in bed, allowing for immediate care. Pads sit on top of the fitted bed sheet, and can also be used on chairs and wheelchairs.

- Collects >1100 ml of fluid
- Superior wicking material helps draw moisture away from patient
- Easy to connect and disconnect
- Available for wired and wireless Parasol Medical™ Monitors
- Integrates with existing nurse call systems

\*Chiarelli PE, Mackenzie LA, Osmotherly PG. Urinary incontinence is associated with an increase in falls: a systematic review. Aust J Physiother. 2009;55(2):89-95. Review. PubMed PMID: 19463079.



## Ultimate Protection by MicrobeCare™

System components are treated with MicrobeCare, a patented, permanently bonded antimicrobial for maximum infection control.

- EPA- and FDA-approved
- Effective against mold, mildew, bacteria and algae
- Provides freshness
- Prevents deterioration and discoloration

## Record a Personalized Message to Patients

A familiar voice can have a powerful effect on a patient who's alone. With the Parasol Wireless Mobility Monitor, family or staff can easily record a personalized message that instructs a patient to sit back down until help arrives. The message will play whenever a patient relieves pressure from any of the pads or unlatches the wheelchair seatbelt.

## Exceptional Versatility

The Parasol Wireless Mobility Monitor can be mounted on a wall, chair or wheelchair. Pressure pads can be placed on any surface where the patient sits or lies. Should a patient get up before your staff can stop them, LED lights allow them to safely navigate from the bed to the bathroom, for additional peace of mind.

## From the Healthcare Innovator

Parasol Medical is a premier developer of specialty medical devices designed to serve the ever-changing healthcare industry. Our greatest priority is making the lives of healthcare workers easier and improving patient outcomes.

At the forefront of new and emerging technologies, Parasol Medical is rapidly expanding its footprint throughout the healthcare market to include products in our four main divisions: Patient Safety, Life Science, Diagnostic and Advanced Wound Care. Parasol Medical complies with Current Good Manufacturing Practices (cGMPs) and maintains an ISO 13485 quality system.

See the Parasol Medical Wireless Patient Fall Prevention System in action at [parasolmed.com/total-room-solution](http://parasolmed.com/total-room-solution).

To learn more about protecting your patients and organization with the Parasol Medical Wireless Patient Fall Prevention System, contact your Parasol Medical sales representative or visit [parasolmed.com](http://parasolmed.com).



“

We knew we had engineered and developed something completely new and innovative when we were awarded 1st place in MedAssets' New Emerging Technology Forum, including more than 500 of the country's most innovative medical device manufacturers.”

**Mike Kilcran**, President of Parasol Medical

## Monitor & Pads

Wireless Monitor & Voice Recorder - Stand Alone	WPADS00
Wireless Monitor & Voice Recorder - Wall Mount	WPADS01
Wireless Monitor & Voice Recorder - Chair Mount	WPADS02
Wireless Chair Pad Sensor 30-Day	WC30
Wireless Bed Pad Sensor 30-Day	WB30
Wireless Chair Pad Sensor 90-Day	WC90
Wireless Bed Pad Sensor 90-Day	WB90
Wireless Floor Pad Sensor 30-Day	WF30
Wireless Potty Pad Sensor 30-Day	WP30
Wireless drySmart® Incontinence Pad	WDS-01
Wireless Seatbelt w/Velcro®	WSB

## Accessories

Wireless Nurse Call Interface	WNCI
Replacement Battery Door Monitor	BCAA
Replacement Battery Door Nurse Call	WNCIBCAA
Wall Mount w/Key	BP02
Monitor Key	KEY01
Parasol AA Battery	PANAA
Wireless User Manual	WMAN-00
Wire Chair Mount 3"	WB01
Wire Chair Mount 6"	WB02
3M Dual Lock	DL01
Wall Mount Lock Key	WLK-01
Mounting Screws for Wall Mount	MTMS-01
drySmart® Wireless Incontinence Quick Connector	QC-02
AC Adapter	WAC-01



1602 Barclay Blvd.  
Buffalo Grove, IL 60089

[parasolmed.com](http://parasolmed.com)  
847.520.3000

FDA Registered Manufacturer | FDA Listed Products | ETL Listed UL Approved | EMC Tested | CE

Distributed by: