

ThermoFlo™ Neo

Now! High performance in a new, smaller size.


High Humidification for patients who need it the most.

- Effective Hygroscopic Condensing Humidifier (HCH) provides protection from respiratory heat loss as well as providing clinically effective humidity.
- Unexcelled low resistance to flow
- Critical in the care of special care babies

Available with or without an offset sampling tube which keeps the adapter weight off of the baby.

#6215; ThermoFlo Neo, 50 per case

#6220; ThermoFlo Neo w/ GS Port,
50 per case



Specifications ¹	Humidity V _T 20/25ml 35 BPM (mg H ₂ O / L air)	Resistance @ 7.5/15 l/min (cm H ₂ O)	Dead Space	Tidal Volume Range	Weight	Housing
ThermoFlo™ Neo	32.5 / 31.0	0.35 / 0.95	2.5 ml	5 – 50 ml	3 gms	Clear

1) Data on File

ARC Medical, Inc.
Anesthesia / Respiratory Care

4296 Cowan Road
Tucker, GA 30084
Order 800-950 ARC1 (2721)

CAUTION: Federal Law restricts this device to sale by or on order of a physician.