


ThermoFlo™ Trach Basic

ThermoFlo™ Trach CRV




ThermoFlo™ Trach CRV

ThermoFlo™ Trach CRV O₂

Humidification: 26mg H₂O/L air, 500 V_T
Resistance: 0.25cm H₂O @ 30 LPM
Dead Space: 12 ml



ThermoFlo™ Trach Basic



Adapter Available for Supplemental O₂

- Now with Cough Relief Valve (CRV)
- ICU / ENT / Home Care
- O₂ Humidified before delivery to patient
- Inspiratory gases are optimally humidified
- High moisture output with low resistance to flow
- Adult & Pediatric Patients
- CRV (Cough Release Valve) may also be used for suctioning
- Humidification at the lowest daily cost

Order Info & Description	Box of 30	8 bx of 30
ThermoFlo™ Trach Basic	#6240	#6240E
ThermoFlo™ Trach CRV w/integral O ₂ Adapter	#6241	#6241E
ThermoFlo™ Trach CRV	#6242	#6242E
O ₂ Adapter for Basic Trach (25/box)	#6258	#6258E

ARC MEDICAL, INC.
Anesthesia / Respiratory Care
4296 Cowan Road
Tucker, GA 30084

www.arcmedical.com
Order Desk: 800-950-ARC1 (2721)
Voice Mail: 800-950-2720