

VeinPress™

The Ultra-Soft Disposable Tourniquet



Place white side on the patient's arm

(place white side up on the patient's arm)



Build up compression force by extending

(build up compression force by extending)



Attach the strap in the velcroed area

(attach the strap in the velcroed area)

The Advantages of VeinPress™

- Easy use elastic strip ensures optimal compression force.
- The ultra-soft fabric provides unmatched comfort for the patient. No more bruising, constriction of the skin, or burns due to hairy arms.
- The Micro-Velcro® closing allows for easy fastening and adjustability.
- Single-use: No infection cross-contamination, cleaning, or disinfection required.
- Cost effective and low weight.

***Servicing the Medical
Community for
Over 20 Years***



VeinPress™

The Ultra-Soft Disposable Tourniquet

Recommended Usage

- Blood donation centers
- Rehabilitation
- Isolation wards
- Pediatrics
- Oncology
- Nursing and retirement homes
- ICU
- Military hospital, MASH
- Emergencies, air rescue, first aid
- Epidemic areas



VeinPress™ Products

Item#	Description	Packaging
MTM430	Pediatric 30mm single-use tourniquet	125/box
MTM435	Adult 35mm single-use tourniquet	100/box

 OVERNIGHT SHIPPING AVAILABLE

 INTERNATIONAL SHIPPING

 MINIMAL PACKAGING

We gladly accept:



For additional information please visit
www.MasonTayler.com



Toll Free US and Canada: 1-888-270-7494
International: 001-716-883-8097
Fax: 716-883-8099
Email: info@masontayler.com
1083 Delaware Ave., Buffalo, NY 14209 USA