

Veinlite

No More Guesswork — No More Probing Making Vein Access As Easy As 1-2-3

Proven Veinlite® Technology and the unique C-Shaped design helps secure the vein, provides traction on the skin and facilitates easy needle insertion.

Clinically Proven Technology in a randomized study of 240 patients at *The Children's Hospital Boston*, Veinlite was found to be significantly better at first attempt in IV placement than standard of care in the Pediatric Emergency Department.

LEDX

- Most Powerful LED colors optimized for vein imaging with increased area of illumination.
 - **Orange** LEDs (24) for superficial veins**
 - **Red** LEDs (8) for deeper veins, dark skin
- Disposable clear plastic cover
- Rechargeable lithium-ion battery with charger
- Weight 83gm (including battery)
- Size 102mm x 65mm x 24mm



LED

- LED colors optimized for vein imaging
 - **Orange** LEDs (12) for superficial veins**
 - **Red** LEDs (12) for deep veins, dark skin
- Disposable clear plastic cover
- Rechargeable lithium-ion battery
- Weight 70gm (3 oz)
- Size 9.5 x 5.5 x 2.4 cm
- One Year Warranty

EMS

- LED colors optimized for vein imaging
 - **Orange** LEDs (12) for superficial veins**
 - **Red** LEDs (4) for deep veins, dark skin
- Disposable clear plastic cover
- Snap-on light shield and Pedi adapter**
- Two AA Batteries (Included)
- Weight 35gm (1.5 oz) without batteries.



PEDI

- LED colors optimized for vein imaging
 - **Orange** LEDs (8) for superficial veins**
 - **Red** LEDs (4) for deep veins, dark skin
- Disposable clear plastic cover
- Small 3 volts Lithium Battery
- Two AA Batteries (Included)
- Weight 40gm (1.5 oz) with battery.

For additional information please call:

Mason Tayler Medical

1083 Delaware Ave., Buffalo, NY 14209 USA
716-883-8097 / 888-270-7494 / Fax: 716-883-8099

Email: masontayler@att.net

Visit our Website: masontayler.com