

## Parasol Advanced Fall Prevention System

The first system that combines fall monitoring with real-time incontinence monitoring.



  
**PARASOL**  
Treated with MicrobeCARE™

# Prevent falls from every surface patients occupy with Parasol Medical's comprehensive solution.

The Parasol Medical Patient Fall Prevention System includes a Bed Pad, Chair Pad, Potty or Commode Pad, Mobility Monitor with Voice Recorder (shown with wall bracket), Nurse Call Interface, drysmart® Incontinence Pads, drysmart™ Quick Connector and a Seat Belt for wheelchair use.



The Parasol Medical Patient Fall Prevention System is the first patient fall prevention aid that combines fall monitoring with real-time incontinence monitoring, for proactive fall prevention.

This comprehensive solution provides caregivers with the most advanced tools to prevent costly injuries due to falls and incontinence. Benefits

include the ability to connect a bed and chair pad to a single monitor. This important safety feature eliminates the potential risk of forgetting to reconnect a pad when transferring a patient between locations.

By Simultaneously monitoring pressure with incontinence and easy integration to the nurse call station-it's the most advanced way to reduce your organization's exposure to risk and protect your patients like never before.



## Revolutionary, Real-Time Incontinence Monitoring

Elimination-related falls are responsible for more than **50%** of falls within acute care facilities.\*

drysmart® incontinence pads monitor wetness in real-time and give nurses more time to reach patients before they get up and risk a fall. Each highly-absorbent pad is engineered with moisture-sensitive thermal conductive ink, and sends an alarm in real-time when a patient evacuates in bed, allowing for immediate care. Pads sit on top of the fitted bed sheet, and can also be used on chairs and wheelchairs.

- Collects >1100 ml of fluid
- Superior breathable/wicking material helps draw moisture away from patient
- Easy to connect and disconnect
- Available for wired and wireless Parasol Medical™ Monitors
- Integrates with existing nurse call systems

\*Chiarelli PE, Mackenzie LA, Osmotherly PG. Urinary incontinence is associated with an increase in falls: a systematic review. Aust J Physiother. 2009;55(2):89-95. Review. PubMed PMID: 19463079.



## Ultimate Protection by MicrobeCare™

System components are treated with MicrobeCare, a patented, permanently bonded antimicrobial for maximum infection control.

- EPA- and FDA-approved
- Effective against mold, mildew, bacteria and algae
- Provides freshness
- Prevents deterioration and discoloration

## Record a Personalized Message to Patients

A familiar voice can have a powerful effect on a patient who's alone. With the Parasol Mobility Monitor, family or staff can easily record a personalized message that instructs a patient to sit back down until help arrives. The message will play whenever a patient relieves pressure from any of the pads or unlatches the wheelchair seatbelt.

## Exceptional Versatility

The Parasol Mobility Monitor can be mounted on a wall, chair or wheelchair. Pressure pads can be placed on any surface where the patient sits or lies. If a patient is incontinent, you can use a drysmart® incontinence pad in conjunction with a pressure pad for additional piece of mind.



Status: Visual Confirmation of active pads.  
Hold: Allows time to move or change a patient.  
Rest: Silence alarm and play back a recorded message.



Long lasting 9V battery easily fits into drop in battery compartment. Select local alarm modes & tones

To learn more about the Parasol Medical Patient Fall Prevention System, contact your Parasol Medical sales representative or visit [parasolmed.com](http://parasolmed.com).



## From the Healthcare Innovator

Parasol Medical is a premier developer of specialty medical devices designed to serve the ever-changing healthcare industry. Our greatest priority is making the lives of healthcare workers easier and improving patient outcomes. At the forefront of new and emerging technologies, Parasol Medical is rapidly expanding its footprint throughout the healthcare market to include products in our four main divisions: Patient Safety, Life Science, Diagnostic and Advanced Wound Care. Parasol Medical complies with Current Good Manufacturing Practices (cGMPs) and maintains an ISO 13485 quality system.

### Monitor & Pads

Monitor & Voice Recorder-Stand Alone	PADS-00
Monitor & Voice Recorder-Wall Mount	PADS-01
Monitor & Voice Recorder-Chair Mount	PADS-02
Bed Pad Sensor 30 Day White 8ft Cord	SB6L8
Bed Pad Sensor 30 Day White 12ft Cord	SB6L12
Bed Pad Sensor 30 Day White 18" Cord	SB6S15
Bed Pad Sensor 90 Day Blue 8ft Cord	SB6L890
Chair Pad Sensor 30 Day White 8ft Cord	SC6L8
Chair Pad Sensor 30 Day White 3ft Cord	SC6S3
Chair Pad Sensor 30 Day White 18" Cord	SC6S15
Chair Pad Sensor 90 Day White 3ft Cord	SC6S390
Potty Pad Sensor 30 Day White 3ft Cord	SP6S3
Wheel Chair Seat Belt	SSBWV

### Accessories

Nurse Call Cord 10ft	25106
Nurse Call Button 10ft Cord	NCB00
Dual Nurse Call Interface	2NCI
Two Pad Y-Splitter	2PI
Replacement Batter Door	BC9V
Wall Mount	BP-01
Locking Wall Mount and Key	BP-02
Lock Key	KEY-01
Chair Mount 3"Width	WB-01
Chair Mount 6"Width	WB-02
Mounting Screws for Wall Mount	MTMS-01
3M Dual Lock for Wall Mount	DL-01
AC Adaptor for Monitor	WAC-01
drysmart® Quick Connect 8ft Cord	QC-03
drysmart® 60x90cm Incontinence Monitoring Pad	WDS-6090
drysmart® 80x90cm Incontinence Monitoring Pad	WDS-8090



Ask about our completely wireless fall prevention system.



1602 Barclay Blvd.  
Buffalo Grove, IL 60089

[parasolmed.com](http://parasolmed.com)  
847.520.3000

FDA Registered Manufacturer | FDA Listed Products | ETL Listed UL Approved | EMC Tested | CE

Distributed by: