

## SternaBra

**SternaBra**, is a post-operative medical Bra, used for maximum sternum wound support after heart surgery

### Features:

- Adjustable shoulder and side straps.
- Breathable woven cotton materials.
- Cross-stitching for steady support.

### Sizes/Ordering codes:

- # SternaBra small
- # SternaBra regular
- # SternaBra large

### Advantages:

- Softer than a normal Bra.
- Steady support to breasts and sternum wound.
- Easily adjustable.
- Front opening for direct access to the wound.
- Offered in three sizes.
- Combinable with SternaSafe Uni.



| Size           | Measure top sternum  | Cup. Size breasts | Order Code |
|----------------|----------------------|-------------------|------------|
| <b>Small</b>   | <b>70cm - 100cm</b>  | <b>A, B</b>       | <b>SBS</b> |
| <b>Regular</b> | <b>90cm - 120cm</b>  | <b>B, C, D</b>    | <b>SBR</b> |
| <b>Large</b>   | <b>110cm - 140cm</b> | <b>D, E, F</b>    | <b>SBL</b> |